



IOL Series Filter Housings

- Single Cartridge Housing Accepts AB Style, Single Open End Filter Elements
- All Stainless Steel Materials of Construction
- Accepts 5" (12.7 cm), 10" (25.4 cm), 20" (50.8 cm) or 30" (76.2 cm) Code 3, Code 7, Code 8 or Code 18 Filter Cartridges
- In-Line Inlet and Outlet Connections for Easy Installation
- Quick Opening T-Handle V-Band Clamp Facilitates Filter Change Out
- Wide Choice of Inlet/Outlet Sizes and O-ring Seal Options
- Passivated for Purity

Housing Specifications

Design Ratings:

260 psig (17.9 bar) @ 100°F (40°C)

200 psig (13.8 bar) @ 300°F (150°C)

NOTE: Above ratings apply to housing only. Cartridge and housing O-ring selection may impose different and narrower limitations.

Materials of Construction:

Head: 316L Stainless Steel

Bowl: 316 Stainless Steel

Connections:

Inlet/Outlet: 3/4" NPT, 1" NPT, 3/4" BSPT, 1" BSPT

Flange Option: 300 lb. ANSI raised face flanges available for 1" IOL Housings

Vent/Drain: 1/4" NPT, BSPT or BSPP


Shell O-rings:

Viton¹ A, FEP Encapsulated Viton A, Silicone Elastomer, Buna N, Ethylene Propylene, Ethylene Propylene for Steam Service




Bracket Option: Stainless steel brackets are available as a standard option.

¹ - Registered trademark of DuPont Dow.




For liquids other than water, multiply differential pressure by specific gravity.



Dimensional Data

Model	Flow ² Rate (gpm/lpm)	Dimensions				Volume Gal/L	Weight Empty	
		A in/mm	B in/mm	C in/mm	D in/mm		Threaded Nozzle lbs/Kg	Flanged Nozzle lbs/Kg
IOL305	5/19	7.8/198	4.81/122	9.81/250	4/102	0.11 /0.4	5.1/2.32	11.9/5.4
IOL31	10/38	16.5/419	4.81/122	9.81/250	12.63/321	0.42/1.6	5.95/2.7	12.8/5.8
IOL32	20/76	26.25/667	4.81/122	9.81/250	22.5/572	0.82/3.1	6.92/3.14	13.7/6.2
IOL33	30/114	36.13/918	4.81/122	9.81/250	32.38/822	1.22/4.6	7.89/3.58	14.7/6.7
IOL705	5/19	7.8/198	4.81/122	9.81/250	4/102	0.11/0.4	5.1/2.32	11.9/5.4
IOL71	10/38	16.5/419	4.81/122	9.81/250	12.63/321	0.42/1.6	5.95/2.7	12.8/5.8
IOL72	20/76	26.25/667	4.81/122	9.81/250	22.5/572	0.82/3.1	6.92/3.14	13.7/6.2
IOL73	30/114	36.13/918	4.81/122	9.81/250	32.38/822	1.22/4.6	7.89/3.58	14.7/6.7
IOL805	5/19	7.8/198	4.81/122	9.81/250	4/102	0.11/0.4	5.1/2.32	11.9/5.4
IOL81	10/38	16.5/419	4.81/122	9.81/250	12.63/321	0.42/1.6	5.95/2.7	12.8/5.8
IOL82	20/76	26.25/667	4.81/122	9.81/250	22.5/572	0.82/3.1	6.92/3.14	13.7/6.2
IOL83	30/114	36.13/918	4.81/122	9.81/250	32.38/822	1.22/4.6	7.89/3.58	14.7/6.7
IOL1805	5/19	7.8/198	4.81/122	9.81/250	4/102	0.11/0.4	5.1/2.32	11.9/5.4
IOL181	10/38	16.5/419	4.81/122	9.81/250	12.63/321	0.42/1.6	5.95/2.7	12.8/5.8
IOL182	20/76	26.25/667	4.81/122	9.81/250	22.5/572	0.82/3.1	6.92/3.14	13.7/6.2
IOL183	30/114	36.13/918	4.81/122	9.81/250	32.38/822	1.22/4.6	7.89/3.58	14.7/6.7

² - Dependent on inlet/outlet selection.



Part Numbers/Ordering Information

IOL ■ ● G ◆ ■ (e.g. IOL 3 1 G N12 H13)

Code ■	Cartridge Style	Code ●	Bowl Length
3	Code 3	05	5"
7	Code 7	1	10"
8	Code 8	2	20"
18	Code 18	3	30"

Code ◆	Inlet/Outlet And Vent/Drain Size
N12	3/4" NPT/1/4" NPT
N16	1" NPT/1/4" NPT
B12	3/4" BSPT/1/4" BSPT
B16	1" BSPT/1/4" BSPT
B17	1" ANSI Flange/1/4" BSPT
N17	1" ANSI Flange/1/4" NPT
P17	1" ANSI Flange/1/4" BSPP
P19	1" BSPP/1/4" BSPP
B35	1" BS4504 Flange/1/4" BSPT
P35	1" BS4504 Flange/1/4" BSPP

Code ■	Seal Materials
H	Viton A
H1	FEP Encapsulated Viton A
H4	Silicone
H13	Buna N
J	Ethylene Propylene
J7	Ethylene Propylene For Steam Service

Visit us on the Web at www.pall.com

Pall Corporation has offices and plants throughout the world in locations including: Argentina, Australia, Austria, Belgium, Brazil, Canada, China, France, Germany, Hong Kong, India, Indonesia, Ireland, Italy, Japan, Korea, Malaysia, Mexico, the Netherlands, New Zealand, Norway, Poland, Puerto Rico, Russia, Singapore, South Africa, Spain, Sweden, Switzerland, Taiwan, Thailand, United Kingdom, United States, and Venezuela. Distributors are located in all major industrial areas of the world.



Pall Corporation

New York - USA
1.800.FILTERS toll free in the USA
516.484.5400 phone
516.484.3216 fax
information@pall.com e-mail

Portsmouth - Europe
44 (0)23 9230 3303 phone
44 (0)23 9230 2510 fax
BioPharmUK@europe.pall.com e-mail

Beijing - China
86 10 87696388 phone
86 10 87691787 fax
pcnhsginfo@pall.com e-mail

Filtration. Separation. Solution.SM

© Copyright 2004, Pall Corporation. Pall, and are trademarks of Pall Corporation. ® Indicates a Pall trademark registered in the USA. *Filtration. Separation. Solution.SM* is a service mark of Pall Corporation.